



Improving consistency for a Mefenamic acid immediate release formulation

Elke Prasad^{a,b,*}, John Robertson^{a, b}, Gavin W. Halbert^{a, b}

^aEPSRC Future Manufacturing Research Hub, University of Strathclyde, Technology and Innovation Centre, 99 George Street, Glasgow, G1 1RD, UK

^bStrathclyde Institute for Pharmacy and Biomedical Sciences, University of Strathclyde, 161 Cathedral Street, Glasgow, G4 0RE, UK

*elke.prasad@strath.ac.uk

Introduction

Mefenamic acid (MFA) product performance is highly dependant on particle size and can lead to variable efficacy.

In this study, a structured development program was applied, developing a crystalline solid dispersion (CSD) formulation, retaining API properties and improving performance consistency compared to a commercial product.




Figure 1: USP II dissolution test (pH 9) of 250mg MFA capsules, Pharmavit Limited (PVL), Batch 4348: red dashed line 85% drug release.

	HUB target product profile
Dose form	Oral solid dose form
Dose	250 mg Mefenamic acid (MFA)
Release profile	<ul style="list-style-type: none"> IR release profile, controlled via particle size Consistent performance
Formulation	Crystalline solid dispersion (CSD)
Process	HME – 3D print/ injection mould
Predictive approach	Rheology screening tool to predict HME process parameters

Table 1: HUB acid target product profile.

Rheology screening




Figure 2: Complex viscosity versus angular frequency for Soluplus +15% w/w Sorbitol (SOL15) at 120°C (blue), 130°C (orange) and 140°C (black).




Figure 3: Storage (filled) and Loss (open) modulus versus angular frequency for SOL15 at 120°C (blue), 130°C (orange) and 140°C (black).




Figure 4: 50% w/w MFA in SOL15SORB after rheology test

Predicted HME process temperature: 132°C

Dominating viscoelastic properties with angular frequency ↑ :

- viscous (G'') >140°C
- elastic (G') <130°C

Screening failed >20% w/w

Hot-melt-extrusion




Figure 5: HME process data. A) Die pressure, B) Torque, C) Rise in product temperature and D) Filament diameter versus screw speed: SOL15: black-140°C, orange-130°C; 50MFA-SOL15: purple-125°C.

Lowest HME process Temperature:

SOL15: 130°C
50MFA-SOL15: 125°C

Product temperature ↑:
SOL15 130°C > 140°C

Die swell SOL15:
130°C > 140°C

Product characterisation




Figure 6: HME extruded filaments: top 50MFA-SOL15, bottom SOL15.




Figure 7: XRPD data of extrudates: top 50MFA-SOL15, bottom SOL15.




Figure 8: SEM image of 50MFA-SOL15 and SOL15 extrudates.





Figure 9: Mechanical properties of 50MFA-SOL15 and SOL15 extrudates: A) Flexural modulus, Strain at max stress, B) Max stress and Modulus of toughness (MoT).





Figure 10: USP II dissolution test (pH 9) of 250mg MFA capsules with 50MFA-SOL15 extrudate. red dashed line 85% drug release.

Conclusion

Rheology guided process development was successful, but also highlights the limitations of this technique.

The targeted immediate release profile for the CSD of MFA in Soluplus®-Sorbitol formulation has been met, whilst the stable crystalline form (I) of MFA has been retained and the consistency of drug release has been significantly improved.

References:

1) Kolter K, Karl M, Gryczke A. 2012. Hot-Melt Extrusion with BASF Pharma Polymers. 2nd Revised and Enlarged Edition ed., Ludwigshafen, Germany: BASF The Chemical Company

Acknowledgements: EPSRC (Grant ref EP/P006965/1), UKRPIF award from HEFCE (Grant ref HH13054). G. W. Halbert is funded by Cancer Research UK (C149/A20496).